


# HOW TO PAIR THE TEE2 PLUS BLOOD GLUCOSE METER WITH THE SMARTLOG APP ON YOUR MOBILE PHONE



## 1. On Phone

Download SmartLog App from either the Apple Store or Google App. Ensure that the Bluetooth setting on your phone is switched on.



## 2. On Meter

Switch on by holding the ● button for 3 seconds until 'SET' appears on the screen.


## 3. On Meter

Press ▼ button once, and then press ● to confirm the change from 'No' to 'Yes' in bottom left of the screen.

## 4. On Meter

The screen will display 'bT'. Press the ▼ button. The meter screen will display 'OFF', 'On', and 'PAIR' in turn. To pair your meter with your smartphone, press the ● button when 'PAIR' blinks on the screen.

### Note:

The  symbol will appear on the meter screen when the Bluetooth feature is on.

## 5. On Phone

Open the SmartLog App.

Press '**More**' at right hand base of screen.



## 6. On Phone

Select '**Accessories**' and then select the TEE2 PLUS meter from the list of devices as the meter to be paired.



## 7. On Phone

Following the steps on the screen;

Press '**Register**'

Select '**Next**'

Select '**Search**'

## 8. On Phone

Select your meter on the SmartLog screen. This is identified by the last 4 digits of the meter serial number which is on the back of the meter.



## 9. On Meter

The meter will now be displaying a six-digit PIN number

## 10. On Phone

Enter this six-digit PIN number into the screen on the SmartLog mobile app and press Pair.

## 11. On Phone

The meter will display 'SUCCESS' once your meter and smartphone are paired and connected.



### Please Note:

If, after pairing, the blood glucose readings are expressed as mg/dl in the SmartLog App on your phone, this can be changed in the '**More**' menu by selecting 'User Information' and changing the Glucose units to mmol/L

